## PLATINUM ENERGY

Platinum Energy Lamp

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Lampu Ajaib Platinum

Patented Ceramic Semiconductors

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Semikonduktor seramik berpaten

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4-14 Microns Energy  
Certified (Japan)

## SAFETY INSTRUCTION

1. The appliance is not to be used by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction.
2. The elderly, children, and those who need medical treatment, must consult with their physicians first before using the appliance.
3. Children being supervised not to play with the appliance.
4. The appliance is not to be used by:
  - persons with heat intolerance
  - persons with poor circulation or peripheral circulation
  - persons under the influence of alcohol
  - persons under doctor's instructions not to use heat lamp
  - infants
  - persons suffering from sunburn
  - persons under medical care for diseases that involve photosensitivity
  - persons receiving photosensitising medications
5. People with diabetes will be more susceptible to burns, please consult a professional before use and use as directed.
6. People with a deep cut or any fresh wounds, please seek immediate medical attention.
7. For wound healing, please use it on the muscles next to the wound, and do not directly use on the wound.
8. If unexpected side effects, such as itching, occur within 48 h of the first session of using an appliance, medical advice should be sought prior to further exposure.
9. In order to avoid overheating, do not cover the heater.
10. Heater is not located immediately below a socket-outlet.
11. Do not use this heater near wet areas in case of electrocution.
12. Do not use this heater in the immediate surroundings of a bath, a shower or a swimming pool.
13. Do not use this heater if it has been dropped.
14. Do not use if there are visible signs of damage to the heater.
15. Use this heater on a horizontal and stable surface.
16. In the process of operation, DO NOT touch the metal screen, as it is hot.
17. WARNING: Do not use this heater in small rooms when they are occupied by persons not capable of leaving the room on their own, unless constant supervision is provided.
18. WARNING: To reduce the risk of fire, keep textiles, curtains, or any other flammable material a minimum distance of 1 m from the air outlet.

### Caution

DO NOT attempt to change or modify any part of the appliance.

DO NOT touch the appliance while the hand is wet.

It is forbidden and dangerous to alter the specification or modify this appliance in any way.

The temperature of the metal screen shaft is between 73°C to 85°C. The distance between the transmitter and the user depends on the individual reaction of acceptable heat. The recommended distance is approximately 25cm from your skin to the Lamp face.

## Guidelines on usage

Wound

Skin

Joint

Stomach

Waist

Head

Eyes

Nose

Throat

Neck & Shoulder

- As with other electrical appliances, please follow the instructions when using.
- Turn on the switch to activate the appliance. After using, turn off the switch immediately so as to save energy and prolong the life span of the appliance.
- The recommended use time is 15 minutes. User should adjust the use time based on individual condition.
- Keep the eyes firmly shut when use on the face.
- Decrease the time per session if the red spots on skin occur.
- If your skin itches, avoid scratching and take a rest.
- If you feel uncomfortable, please stop using the appliance and consult with doctor immediately.

- The appliance has a protective device for overheating. The electric current will be automatically cut-off once overheated. After cooling, the electric current will resume automatically.
- There should be an accompanying adult on hand when an elderly person/chronic illness sufferer desires to use the appliance.
- User with special medical problem should consult the doctor before using.
- Patients with cardiovascular and chronic diseases must use the appliance only under the supervision of a competent adult.

In case of any abnormal situation, please switch off the appliance; remove the plug and contact your distributor.

Follow this user manual book strictly and do not use any parts outside the specifications.

## Guidelines on maintenance

- Beware of the following safety precautions to avoid accident
  - \* When the appliance is in operation, do not twist, pull or blind the electric cord. Do not put any heavy object on the cord.
  - \* Use only outlet of AC 220-240V
  - \* Once the cord torn or faulty, do not use the appliance and contact the distributor for repair.
- For cleaning and maintenance, use soft and damp cloth or mild detergent to wipe. DO NOT use volatile cleaning liquid.
- After a period of usage, use a hair dryer (use only cool wind) to spout out the dirt in the transmitter.

### Platinum Energy Lamp Specifications

Power Requirement	: AC220-240V / 50Hz
Power Output	: 200W±10%
Size	: 14cm(W) X 38.5cm(H)
Weight	: approx. 1.65kg

This product conforms CE safety standard.

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